

Pedestrian

Model Number: PED



“Sit and Stand” Pedestrian versions of the popular Hybrid-II 50th, Hybrid-III 50th Male, 95th Large Male, and the 5th Small Female are available from DENTON ATD, INC. Because of their unique design, Pedestrian dummies have a varied range of applications including assessment of vehicle component interaction and vehicle contact points. These dummies are complete and certified with all standard upper torso calibration (according to model) allowing for head and thorax injury assessment. The pedestrian series of dummies are equipped with a straight lumbar spine, a modified pelvic structure and knee slider replacements.

Hybrid III Pedestrian **Includes:**

- ◆ Chest Deflection Transducer
- ◆ Six-Axis Neck Transducer Structural Replacement
- ◆ Accelerometer Mountings for Head, Chest, and Pelvis. User specified at time of order placement

**Other instrumentation available but not included

Hybrid II Pedestrian **Includes:**

- ◆ Accelerometer Mounts for the Head and Chest. User specified at time of order placement

Hybrid III 50th Pedestrian Components.

78051-60P	Pelvis (Pedestrian)
78051-41P	Hip Flesh, Left
78051-42P	Hip Flesh, Right
78051-53P	Lumbar Spine Bracket
78051-88P	Lumbar to Thoracic Adaptor
78051-66P	Lumbar Spine
78051-24P	Knee Slider Replacement, Inboard
78051-25P	Knee Slider Replacement, Outboard

SA-150-M002PED Hybrid II 50th Percentile Male Pedestrian Test Dummy

78051-218-PED Hybrid III 50th Percentile Male Pedestrian Test Dummy

880995-000PED Hybrid III 95th Percentile Male Pedestrian Test Dummy

880105-1PED Hybrid III 5th Percentile Female Pedestrian Test Dummy

For more information about this test device or test dummy configurations specific to your needs, please contact the Customer Service Department at DENTON ATD or log onto www.dentonatd.com.